



WWW.DENTALLEARNING.NET

DENTAL LEARNING

A PEER-REVIEWED PUBLICATION

Knowledge for Clinical Practice

OPTIMIZING LIGHT-CURING OUTCOMES IN RESTORATIVE DENTISTRY

By Scott Coleman, DDS, MAGD

INSIDE
Earn 2
CE
Credit

Written for
dentists, hygienists,
and assistants



Approved PACE Program Provider
FAGD/MAGD Credit Approval does not imply acceptance by a
state or provincial board of dentistry or AGD endorsement.
2/1/2024 - 1/31/2030
Provider ID: # 346890
AGD Subject Code: 250



Dental Learning, LLC is a Dental Board of California CE Provider.
The California Provider # is RP5062. All of the information contained on this certificate is truthful and
accurate. Completion of this course does not constitute authorization for the attendee to perform any
services that he or she is not legally authorized to perform based on his or her license or permit type.
This course meets the Dental Board of California's requirements for 2 units of continuing education.
CA course code is 01-5062-25316



ABSTRACT

This continuing dental education print supplement provides a comprehensive review of light-curing technology principles, clinical best practices, and evidence-based considerations for optimizing composite resin polymerization in modern restorative dentistry. Using the Kerr Demi™ Pro Curing Light as a clinical reference platform, participants will explore the science of photopolymerization, the significance of irradiance output and beam collimation, ergonomic factors influencing clinician fatigue and procedural accuracy, and the role of integrated safety technologies in producing consistent, high-quality restorations.

LEARNING OBJECTIVES

Upon completion of this course, participants will be able to:

- Describe the fundamental principles of photopolymerization and the role of light irradiance in achieving adequate composite cure depth and hardness.
- Identify the clinical implications of inadequate light-curing, including marginal failure, secondary caries, and restoration longevity concerns.
- Explain how collimated beam output and curing tip design affect irradiance delivery at clinical distances.
- Evaluate the safety benefits of Periodic Level Shift (PLS) technology in preventing overheating and protecting pulpal and soft tissue integrity.
- Discuss ergonomic design considerations in curing light selection and their relationship to clinician fatigue and procedural efficiency.
- Apply knowledge of multimode curing protocols (Standard, High, Ortho) to appropriate clinical scenarios.
- Recognize the diagnostic applications of transillumination in early-caries detection and fracture assessment.

ABOUT THE AUTHOR



Scott Coleman, DDS, MAGD

Dr. Coleman graduated cum laude from Stephen F. Austin State University in 1980 with a major in biology and a minor in chemistry. He earned his Doctor of Dental Surgery degree in 1984, graduating first in his class from the University of Texas Dental Branch at Houston. A member of the AGD, AACD, ADA, and Texas Dental Association, Dr. Coleman has attained Fellowship standing in the AGD, International Academy of Dental Facial Esthetics, American College of Dentists, International College of Dentists, and the International Congress of Oral Implantologists. He is currently a clinical Adjunct Professor with the University of Texas Health Science Center at Houston School of Dentistry and serves on the Board of Regents for his alma mater, Stephen F. Austin State University.

SPONSOR/PROVIDER: This is a Dental Learning, LLC continuing education activity. **STATEMENTS:** Dental Learning, LLC designates this activity for 2 CE credits. Dental Learning, LLC is also designated as an Approved PACE Program Provider by the Academy of General Dentistry. The formal continuing education programs of this program provider are accepted by AGD for Fellowship, Mastership, and membership maintenance credit. Approval does not imply acceptance by a state or provincial board of dentistry or AGD endorsement. The current term of approval extends from Feb 1, 2024 to Jan 31, 2030. Provider ID: # 346890. **EDUCATIONAL METHODS:** This course is a self-instructional journal and web activity. Information shared in this course is based on current information and evidence. **PUBLICATION DATE:** June 2026. **EXPIRATION DATE:** June 2029. **REQUIREMENTS FOR SUCCESSFUL COMPLETION:** To obtain 2 CE credits for this educational activity, participants review the material, complete the course evaluation and obtain a score of at least 70%. **AUTHENTICITY STATEMENT:** The images in this course have not been altered. **SCIENTIFIC INTEGRITY STATEMENT:** Information shared in this continuing education activity is developed from clinical research and represents the most current information available from evidence-based dentistry. **KNOWN BENEFITS AND LIMITATIONS:** Information in this continuing education activity is derived from data and information obtained from the reference section. **EDUCATIONAL DISCLAIMER:** Completing a single continuing education course does not provide enough information to result in the participant being an expert in the field related to the course topic. It is a combination of many educational courses and clinical experience that allows the participant to develop skills and expertise. **PROVIDER DISCLOSURE:** Dental Learning, LLC does not have a leadership position or a commercial interest in any products that are mentioned in this article. No manufacturer or third party has had any input into the development of course content. **CE PLANNER DISCLOSURE:** The planner of this course, Jaci Soukhathivong, does not have a leadership or commercial interest in any product or services discussed in this educational activity. Any questions or comments can be sent to Support@dentallearning.net. **TARGET AUDIENCE:** This course was written for dentists, dental hygienists, and assistants, from novice to skilled. Please direct all questions pertaining to Dental Learning, LLC or the administration of this course to Support@dentallearning.net. Go Green, Go Online to <https://www.dentallearning.net/article/light-curing-outcomes> to take this course. © 2026



Section 1:

The Science of Photopolymerization in Restorative Dentistry

Light-cured composite resins have become the restorative material of choice in contemporary dentistry due to their excellent esthetic properties, mechanical performance, and clinical versatility. The ability to control the setting reaction through light exposure gives clinicians command over working time and polymerization kinetics, a significant clinical advantage over chemically cured materials. However, the quality of polymerization is not guaranteed by material chemistry alone—it is fundamentally dependent on the quality of the light delivered to the restoration.

Photopolymerization is initiated when light of an appropriate wavelength activates photoinitiators within the resin matrix. The most commonly used photoinitiator in contemporary composites is camphorquinone (CQ), which absorbs light in the blue wavelength range of approximately 400–500 nm, with peak absorption near 468 nm. Upon light absorption, CQ generates free radicals that initiate chain polymerization of the methacrylate monomers. More recent formulations also incorporate alternative photoinitiators such as phenyl-propanedione (PPD) and acylphosphine oxides (APO), which may require broader or slightly shifted wavelength coverage for optimal activation.

The degree of conversion (DC), the proportion of carbon-carbon double bonds converted from monomer to polymer, is the primary determinant of restoration quality. Higher degrees of conversion correlate directly with superior mechanical properties, including flexural strength, hardness, and wear resistance, as well as reduced monomer elution and improved biocompatibility. Incomplete polymerization can leave residual monomer in the restoration matrix, which may contribute to cytotoxicity, microleakage, secondary caries, and premature restoration failure.

Irradiance: The Critical Variable

The single most important physical parameter governing the degree of conversion is irradiance, measured in milliwatts per square centimeter (mW/cm^2). Irradiance

is the radiant power delivered per unit area, and it determines the rate and completeness of photoinitiator activation. Clinical studies consistently demonstrate that higher irradiance values produce superior polymerization outcomes, particularly in the deeper layers of thick restorations, through opaque shades, or when curing at increased working distances from the tooth surface.

The American Dental Association and numerous clinical researchers recommend a minimum irradiance of $400 \text{ mW}/\text{cm}^2$ for effective composite cure, with higher values, often exceeding $1,000 \text{ mW}/\text{cm}^2$, preferred for fast-cure protocols and optimal material performance. Critically, irradiance declines with increasing distance between the light tip and the restoration surface. In many clinical scenarios, direct contact placement of the curing tip is not feasible due to interproximal access limitations, matrix band presence, or anatomical constraints. A curing light that maintains high irradiance over a range of clinical distances is therefore a meaningful clinical asset.

Section 2:

Clinical Consequences of Inadequate Light-Curing

The consequences of under-cured composite restorations extend well beyond the immediate procedural outcome. Inadequate polymerization has been linked to a cascade of clinical complications that compromise patient outcomes, increase the burden of retreatment, and affect the long-term success rates of tooth-colored restorations.

Mechanical and Structural Failures

Incompletely polymerized composite resins exhibit reduced hardness, lower flexural strength, and diminished wear resistance compared to fully cured material. These mechanical deficiencies predispose restorations to surface wear, fracture under occlusal loading, and marginal degradation over time. Clinically, the result is premature restoration failure, replacement cycles that increase patient cost, tooth structure loss, and clinical chair time.

Marginal Integrity and Secondary Caries

Polymerization shrinkage stress at the tooth-restoration interface is a well-documented challenge in



composite dentistry. When resins are inadequately cured, the shrinkage stress is not fully developed, which may paradoxically reduce immediate marginal gap formation. However, the restoration's inferior mechanical properties mean it is more susceptible to microcracking, marginal breakdown, and eventual microleakage over time. Microleakage enables ingress of oral fluids and cariogenic bacteria, a leading cause of secondary (recurrent) caries beneath restorations, one of the most common reasons for restoration replacement in dental practice.

Biocompatibility Concerns

Residual unpolymerized monomers, particularly bisphenol A-glycidyl methacrylate (Bis-GMA), triethylene glycol dimethacrylate (TEGDMA), and hydroxyethyl methacrylate (HEMA), are known to elute from incompletely cured restorations into the oral environment. These monomers have demonstrated cytotoxic and genotoxic potential in laboratory studies, with implications for pulpal biocompatibility and potential sensitization reactions. Maximizing degree of conversion through adequate light-curing is an important biocompatibility strategy.

Given these clinical ramifications, the selection, maintenance, and proper use of a high-performance curing light is not a matter of convenience, it is a core component of evidence-based restorative practice.

Section 3: Engineering the Ideal Curing Light: Beam Collimation and Output Consistency

Modern curing lights vary considerably in their optical design, output characteristics, and clinical performance. Understanding the engineering principles behind a high-performance device allows clinicians to make informed purchasing decisions and to use their equipment in ways that maximize patient benefit.

Collimated vs. Divergent Light Output

Light emitted from an LED curing light can be either collimated (parallel beam) or divergent (spreading outward from the tip). Collimated output maintains intensity over a broader range of distances from the tip, whereas divergent output—while delivering high intensity at point-blank

range, drops off more rapidly as distance increases. This distinction has critical clinical implications.

Consider the reality of posterior composite restorations: A clinician curing a Class II restoration through a matrix band, or reaching a distal surface of a second molar, may be operating at 4–10 mm of effective working distance from the restoration surface. If the curing light delivers collimated output, the irradiance at 6 mm may be 70–80% of its tip-contact value. With a highly divergent beam, irradiance at the same distance may fall to 40–50%, potentially below the threshold for adequate polymerization.

Clinical Data: The Demi Pro delivers 1,393 mW/cm² at direct contact and maintains 1,025 mW/cm² at 6-mm working distance through a tooth, representing a 26% reduction. Competitive devices in the same class show 27% and 39% reductions, respectively, at the same distance (data on file, Kerr Corporation).

This sustained output at clinical distances is the functional signature of collimated beam technology. It translates directly into more reliable curing results in posterior regions, in deep boxes, and in situations where direct tip placement is anatomically or procedurally restricted.



Aperture Size and Beam Profile

The aperture diameter of the curing tip determines the effective curing area. Wider apertures (8–10 mm) cover a larger surface area per exposure but distribute energy over a greater zone, potentially reducing point irradiance. Narrower apertures concentrate light over a smaller area with higher local intensity. Clinical research suggests that a range of aperture sizes is valuable, and the Demi Pro's performance data across 4-mm, 6-mm, and 10-mm apertures demonstrates consistent output curves, ensuring predictable results regardless of which tip configuration

Optimizing Light-Curing Outcomes in Restorative Dentistry

is selected. The 360° rotatable, low-profile curing head design facilitates optimal tip-to-tooth proximity even in anatomically challenging regions such as the distal of second molars, furcation areas in periodontal procedures, and pediatric dentistry applications.

Section 4: Thermal Safety in Light-Curing: PLS Technology and Pulpal Protection

One of the most significant patient safety considerations in operative dentistry is the management of heat generated during light-curing procedures. Elevated temperatures at the pulpal wall of a preparation can induce irreversible pulpitis, pulpal necrosis, or accelerated dentin hypersensitivity, outcomes that are avoidable with appropriate technology and technique.

The Thermal Challenge of High-Output LEDs

High-irradiance curing lights are inherently more powerful heat generators than their predecessors. The same photonic energy that drives rapid polymerization also converts partially to thermal energy at the restoration surface and within the handpiece itself. Conventional high-output LED lights may generate surface temperatures exceeding safe thresholds during prolonged exposures, particularly when curing through opaque materials, applying multiple increments in succession, or working in proximity to thin dentin walls.

Research by various authors has identified a critical temperature rise threshold of 5.6°C at the pulpal wall, beyond which irreversible pulpal damage may occur. Managing heat generation during curing is therefore not merely an equipment engineering concern, it is a patient safety imperative.


Periodic Level Shift (PLS) Technology

The Demi Pro employs a patented thermal management strategy called Periodic Level Shifting (PLS), which represents a meaningful advance in curing light safety engineering. Rather than delivering light at a constant maximum output throughout the curing cycle, which would maximize heat generation, PLS technology modulates the output intensity periodically throughout the exposure.

During a standard 5-second curing cycle, PLS technology shifts the output intensity multiple times, alternating between high-output pulses and brief reduced-output intervals. This pulsed modulation strategy maintains effective radiant energy delivery for polymerization while interrupting the continuous heat accumulation that characterizes constant-output devices. The result is sustained polymerization performance without the thermal load associated with constant-maximum operation.

The PLS output profile during a 5-second exposure shows cyclic shifts between approximately 1,100 and 1,400 mW/cm², delivering an effective curing dose while preventing sustained high-temperature output. This thermal management approach is particularly valuable in situations involving thin dentin walls, vital teeth with recent trauma, direct pulp caps, or patients with pre-existing sensitivity.

The PLS output profile during a 5-second exposure shows cyclic shifts between approximately 1,100 and 1,400 mW/cm², delivering an effective curing dose while preventing sustained high-temperature output. This thermal management approach is particularly valuable in situations involving thin dentin walls, vital teeth with recent trauma, direct pulp caps, or patients with pre-existing sensitivity.





Integrated Radiometer

The Demi Pro includes a built-in radiometer integrated into the charging base, allowing clinicians to verify curing light output before each clinical session. This feature addresses a well-documented clinical problem: Curing lights in active practice frequently deliver substantially less output than their specifications indicate due to LED degradation, tip contamination, or mechanical damage. Studies in dental practice settings have found that a significant proportion of curing lights in clinical use produce irradiance values below recommended minimums, often without the clinician's awareness.

Regular output verification using an integrated radiometer allows for early identification of device degradation, proactive maintenance, and confident procedural planning. A clinician who knows their curing light is performing at specification can make informed decisions about exposure times and incremental curing protocols with greater confidence.



Section 5:

Ergonomics, Clinician Fatigue, and Procedural Efficiency

The ergonomic design of dental instruments has received growing attention in recent years as the profession recognizes the significant burden of musculoskeletal disorders (MSDs) among dental practitioners. Studies indicate that a substantial proportion of dentists and dental hygienists experience work-related musculoskeletal symptoms during their careers, with the shoulder, neck, and wrist being particularly vulnerable areas. Instrument

design plays a meaningful role in reducing ergonomic strain and supporting long-term career sustainability.



Weight and Balance in Curing Light Design

A curing light used dozens of times per day must be genuinely comfortable in hand—not merely acceptable. The Demi Pro was engineered to weigh just 115 grams, placing it among the lightest professional-grade curing lights available. This pen-like weight distribution reduces the sustained static loading on the hand, wrist, and forearm that occurs when clinicians hold the device at extension for multiple curing cycles in succession.

The significance of this weight reduction is best appreciated in a high-volume restorative practice, where cumulative exposure to even modest ergonomic stressors compounds into meaningful physical strain over months and years of practice. Ergonomic instruments that reduce per-procedure stress contribute to practitioner longevity and help maintain the precision of hand positioning through the end of a demanding clinical day.

360° Rotating Tip and Low-Profile Access

The rotating curing head is one of the Demi Pro's most clinically significant design features. Traditional fixed-tip curing lights require the clinician to reposition their wrist and sometimes their entire arm posture to achieve adequate light angulation on posterior, lingual, or distal surfaces. This forced awkward posturing is a known ergonomic risk factor for shoulder and cervical strain.

The 360° rotating low-profile tip allows the clinician to maintain a neutral wrist and arm position while directing

Optimizing Light-Curing Outcomes in Restorative Dentistry



light output to virtually any intraoral surface. The tip can be rotated to the optimal angle before the curing cycle begins, eliminating mid-procedure repositioning and enabling consistent, controlled light delivery without compromising clinician posture. Clinical feedback from practitioners confirms that this feature makes curing at every angle genuinely comfortable, supporting extended clinical use without progressive fatigue.

Curing Modes and Time Efficiency

The Demi Pro offers three clinical curing modes: Standard, High, and Ortho, each with selectable curing times, providing the flexibility to match the curing protocol to the specific material, restoration type, and clinical context. The mode selection is displayed on the device's integrated screen, enabling at-a-glance confirmation of the active protocol.

In high-volume practice settings, the ability to confidently complete a curing cycle in 5 seconds while trusting that the output is sufficient for deep, reliable polymerization translates into real scheduling efficiency. Cumulative time savings across dozens of restorations per day represent a meaningful contribution to practice productivity without compromising quality of care.

Section 6: Transillumination as a Diagnostic Adjunct

Beyond its primary function as a curing device, the

Demi Pro platform supports an optional Transillumination Tip accessory that extends the device's clinical utility into the realm of diagnostic imaging. Transillumination is a noninvasive, radiation-free diagnostic modality that uses transmitted light to visualize internal tooth structure and identify anomalies that may not be detectable by conventional visual examination or bitewing radiography alone.

Principles and Clinical Applications



When intense light is directed through tooth structure, variations in optical density create visible contrast. Healthy enamel transmits light relatively uniformly, while areas of demineralization, carious lesion formation, or structural disruption absorb and scatter light differently, producing detectable shadows or bright zones in the transilluminated image.

The clinical applications of transillumination are broad and clinically meaningful. Interproximal caries detection using fiber-optic transillumination (FOTI) has been validated in peer-reviewed research as a reliable adjunct to radiographic examination, with sensitivity for enamel and dentinal caries that is competitive with conventional bitewing radiography in certain anatomical regions. The technique is particularly valuable in patients for whom radiographic exposure is to be minimized, such as pregnant patients, pediatric patients with high cumulative radiation histories, or patients with documented radiation sensitivity.


Transillumination is also valuable in the assessment of vertical root fractures, craze lines, and occlusal decalcification. Vertical root fractures are among the most diagnostically challenging conditions in clinical dentistry, often producing nonspecific symptoms and limited radiographic findings until late-stage pathology has developed. Transillumination can illuminate fracture lines in ways that guide clinical decision-making and help clinicians counsel patients more accurately about prognosis and treatment options.

Integration with the Demi Pro Platform

The availability of the Transillumination Tip as an accessory for the Demi Pro handpiece creates a clinically efficient platform in which a single device serves both therapeutic and diagnostic functions. The tip accepts the



Interproximal caries



Diagnosis of a vertical fracture and occlusal decalcification



Dianosis of a vertical crack on a tooth undergoing endodontics

same handpiece without modification, and its use requires no additional operator training beyond familiarity with standard transillumination technique.

Section 7: Clinical Protocol Recommendations and Best Practices

Evidence-based use of a high-performance curing light requires attention not only to device selection and settings, but also to the broader clinical protocol context in which light-curing occurs. The following recommendations represent best practices synthesized from clinical research, materials science evidence, and the specific capabilities of the Demi Pro system.

Incremental Layering Technique

Regardless of curing light output, composite restorations should be placed in increments of no greater than 2-mm depth to ensure adequate light penetration and uniform polymerization throughout the increment. Bulk-fill composites specifically formulated for deeper increments may allow up to 4–5 mm placements per increment, but manufacturers' specifications for curing time and minimum irradiance should be strictly followed.

The Demi Pro's collimated output provides a meaningful advantage for consistent incremental curing, ensuring that light intensity at the base of each increment meets polymerization thresholds.

Curing Distance Management

The curing tip should be positioned as close to the restoration surface as clinically feasible, ideally within 1–2 mm. The 360° rotating tip facilitates this by allowing angulation optimization without forced wrist positioning. When a matrix band or anatomical obstruction prevents direct placement, the Demi Pro's performance at extended distances provides a meaningful safety margin.

Output Verification

Clinicians should establish a routine protocol for verifying curing light output using the integrated radiometer in the charging base. Output verification should occur at the beginning of each clinical day and following any period of extended storage, device impact, or replacement of the curing tip. Documenting output values as part of equipment maintenance records supports quality assurance in the clinical setting.

Light Shield Use

The curing light package includes light shields to protect patients and clinical staff from direct blue light exposure during curing procedures. Operator eye protection should be used consistently—the high irradiance values that make the Demi Pro an effective curing instrument also necessitate appropriate protective protocols.

Battery Management

The Demi Pro is supported by multiple long-lasting batteries, with inductive charging through the integrated base. Consistent return of the device to its charging base between cases ensures the battery is maintained at an appropriate charge level. A battery charge indicator on the device display allows clinicians to monitor power status proactively, preventing unexpected output reduction mid-procedure.

Section 8: Comparative Clinical Performance and Evi-

dence Summary

A comprehensive clinical evaluation of curing light performance requires consideration of multiple parameters: irradiance at clinical distances, thermal output, ergonomic characteristics, mode versatility, and longevity of output over the device's service life. Synthesizing these parameters provides a basis for evidence-informed device selection.

Parameter	Demi Pro	VALO™ X	SmartLite® Pro
Mode	Standard	Std Power (Steady Blue)	Standard
Irradiance at 0 mm	1,393 mW/cm ²	1,137 mW/cm ²	1,233 mW/cm ²
Irradiance at 6 mm	1,025 mW/cm ²	830 mW/cm ²	748 mW/cm ²
Thermal Protection	PLS Technology	Standard	Standard
Tip Rotation	360°	Fixed	Fixed
Built-in Radiometer	Yes	No	No

The comparative data demonstrates that the Demi Pro achieves the highest irradiance values at both point-blank and clinical working distances among comparable devices. The maintenance of over 1,000 mW/cm² at 6 mm is particularly notable, as this is a clinically relevant distance for many posterior restoration scenarios. The combination of superior output consistency, thermal management

through PLS technology, ergonomic lightweight design, and integrated diagnostic radiometer positions the Demi Pro as a comprehensive clinical solution.

Conclusion

Light-curing technology is foundational to the performance of contemporary composite resin restorations. The quality of polymerization achieved during a restorative procedure directly determines the mechanical properties, longevity, and biocompatibility of the restoration, making the curing light one of the most consequential instruments in the restorative armamentarium.

The Kerr Demi™ Pro Curing Light represents a thoughtful integration of engineering advances that address the clinical challenges most commonly encountered in light-curing practice: irradiance maintenance at clinical distances, thermal safety management, ergonomic design for reduced clinician fatigue, mode versatility for diverse procedural applications, and diagnostic expansion through transillumination capability.

Clinicians who understand the principles underlying these technologies are better equipped to use them effectively, to make informed comparisons when selecting equipment, and to communicate with patients about the relationship between procedural quality and restoration outcomes. Continuing education in this domain is an investment not only in clinical knowledge, but also in the sustained quality and safety of patient care.

References

- Rueggeberg FA. State-of-the-art: Dental photocuring, a review. *Dent Mater.* 2011;27(1):39–52.
- Leprince JG, Palin WM, Hadis MA, Devaux J, Leloup G. Progress in dimethacrylate-based dental composite technology and curing efficiency. *Dent Mater.* 2013;29(2):139–156.
- Price RB, Shortall AC, Palin WM. Contemporary issues in light curing. *Oper Dent.* 2014;39(1):4–14.
- Daronch M, Rueggeberg FA, Hall G, De Goes MF. Effect of composite temperature on in vitro intrapulpal temperature rise. *Dent Mater.* 2007;23(10):1283–1288.
- Zach L, Cohen G. Pulp response to externally applied heat. *Oral Surg Oral Med Oral Pathol.* 1965;19:515–530.
- Strassler HE, Price RB. A review of FOTI and its clinical application in caries detection. *Compend Contin Educ Dent.* 2014;35(5):344–348.
- Polín V, Vocánková I, Konýšková R, et al. Composite degree of conversion: Influence of light curing mode and distance. *J Dent.* 2020;95:103318.
- Bruzell EM, Johnsen B, Aalerud TN, Dahl JE, Christensen T. In vitro efficacy and risk for adverse effects of light-assisted tooth bleaching. *Photochem Photobiol Sci.* 2009;8(3):377–385.
- Kerr Corporation. Demi™ Pro Curing Light Product Brochure. Data on File. MKT-25-1336 Rev 1. ©2026.
- Tauböck TT, Tarle Z, Marovic D, Attin T. Pre-heating of high-viscosity bulk-fill resin composites: Effects on shrinkage force and monomer conversion. *J Dent.* 2015;43(11):1358–1364.



CE Quiz

To complete this quiz online and immediately download your CE verification document, visit <https://www.dentallearning.net/article/light-curing-outcomes>, then log in to your account (or register to create an account). Upon completion and passing of the exam, you can immediately download your CE verification document.

Optimizing Light-Curing Outcomes in Restorative Dentistry:

A Clinical and Technical Review of the Demi™ Pro Curing Light

Instructions: Select the single best answer for each question. A score of 70% or higher (21 or more correct answers) is required for CE credit. Circle or record your answer on a separate answer sheet before checking the answer key

Section 1: Photopolymerization Science (Questions 1–5)

1. What is the primary photoinitiator used in most contemporary composite resin formulations?
 - A. Phenylpropanedione (PPD)
 - B. Camphorquinone (CQ)
 - C. Acylphosphine oxide (APO)
 - D. Bisphenol A-glycidyl methacrylate (Bis-GMA)
2. The degree of conversion (DC) in a composite restoration refers to:
 - A. The percentage of filler content converted to polymer
 - B. The proportion of carbon-carbon double bonds converted from monomer to polymer
 - C. The ratio of light output to heat generated during curing
 - D. The change in color of the composite from uncured to cured state
3. Which of the following wavelength ranges activates camphorquinone most effectively?
 - A. 350–390 nm (UV range)
 - B. 520–560 nm (green range)
 - C. 400–500 nm (blue range)
 - D. 600–650 nm (red-orange range)
4. Irradiance in the context of dental curing lights is measured in:
 - A. Joules per cubic centimeter (J/cm^3)
 - B. Milliwatts per square centimeter (mW/cm^2)
 - C. Nanometers (nm)
 - D. Lux (lx)
5. What minimum irradiance value is generally recommended by clinical researchers for effective composite polymerization?
 - A. 100 mW/cm^2
 - B. 200 mW/cm^2
 - C. 400 mW/cm^2
 - D. 800 mW/cm^2

Section 2: Clinical Consequences of Inadequate Curing (Questions 6–10)

6. Which of the following is a direct mechanical consequence of incomplete composite polymerization?
 - A. Increased surface hardness
 - B. Enhanced flexural strength
 - C. Reduced wear resistance
 - D. Improved marginal seal
7. Secondary (recurrent) caries beneath composite restorations is most directly associated with:
 - A. Excess monomer elution causing pulpal inflammation
 - B. Microleakage enabling ingress of oral fluids and cariogenic bacteria
 - C. Elevated surface temperature during polymerization
 - D. Photoinitiator residue at the tooth-restoration interface

8. Which monomer commonly found in incompletely cured composites has demonstrated cytotoxic and genotoxic potential?
 - A. Camphorquinone
 - B. Hydroxyapatite
 - C. TEGDMA (triethylene glycol dimethacrylate)
 - D. Zinc oxide
9. Polymerization shrinkage stress at the tooth-restoration interface is best managed by:
 - A. Increasing the irradiance to the highest possible setting
 - B. Using bulk increments of 6 mm or greater
 - C. Incremental placement and adequate light-curing of each layer
 - D. Reducing curing time to minimize heat generation
10. A clinician notices a recently placed Class II composite restoration showing marginal staining and early secondary caries 18 months post-placement. Which factor is MOST likely to have contributed to this outcome?
 - A. Shade selection that was too opaque
 - B. Inadequate light-curing leading to microleakage
 - C. Excess finishing and polishing of the restoration margins
 - D. Use of a light-cured vs. dual-cured material

Section 3: Beam Collimation and Irradiance (Questions 11–15)

11. A collimated light beam from a curing light is best described as one that:
 - A. Produces maximum heat output at the tip
 - B. Spreads outward rapidly with increasing distance
 - C. Maintains parallel output with sustained intensity over distance
 - D. Operates exclusively at a fixed wavelength of 468 nm
12. Why is maintaining irradiance at extended clinical distances (4–10 mm) particularly important in posterior restorative dentistry?
 - A. Posterior composites contain different photoinitiators requiring lower irradiance
 - B. Matrix bands, anatomical access, and restoration geometry may prevent direct tip contact
 - C. Blue light is absorbed more rapidly in posterior regions of the oral cavity
 - D. Extended distances are required to prevent thermal damage in molar teeth
13. According to the comparative data presented in this white paper, the Demi Pro's irradiance at 6 mm through a tooth is approximately:
 - A. 748 mW/cm^2
 - B. 830 mW/cm^2
 - C. 1,025 mW/cm^2
 - D. 1,393 mW/cm^2

Optimizing Light-Curing Outcomes in Restorative Dentistry

14. The Demi Pro's irradiance at direct tip contact (0 mm) in Standard mode is reported as:

- A. 748 mW/cm²
- B. 1,137 mW/cm²
- C. 1,233 mW/cm²
- D. 1,393 mW/cm²

15. A wider curing tip aperture (e.g., 10 mm vs. 4 mm) primarily affects light delivery by:

- A. Increasing peak irradiance at the center of the beam
- B. Covering a larger surface area per curing exposure
- C. Producing a more collimated output over distance
- D. Reducing the wavelength emitted by the LED

Section 4: Thermal Safety and PLS Technology (Questions 16–20)

16. The critical temperature rise threshold at the pulpal wall beyond which irreversible pulpal damage may occur is approximately:

- A. 2.0°C
- B. 3.5°C
- C. 5.6°C
- D. 8.0°C

17. PLS (Periodic Level Shift) technology in the Demi Pro manages thermal output by:

- A. Reducing total irradiance to below 800 mW/cm² throughout the cycle
- B. Periodically modulating output intensity during the curing cycle
- C. Using a cooling fan integrated into the handpiece body
- D. Limiting each curing exposure to a maximum of 2 seconds

18. In which clinical scenario would PLS thermal management be MOST critical?

- A. Curing a composite veneer on a fully crowned tooth
- B. Light-curing a pit-and-fissure sealant on a primary molar
- C. Curing a deep restoration with thin residual dentin overlying the pulp
- D. Polymerizing a composite core buildup on an endodontically treated tooth

19. The integrated radiometer in the Demi Pro system is designed to:

- A. Measure the patient's pulpal temperature during curing
- B. Verify the curing light's output before clinical use
- C. Detect residual monomer in the cured composite surface
- D. Calibrate the wavelength spectrum of the LED emitter

20. Studies of curing lights in active dental practice settings have found that:

- A. Nearly all curing lights maintain output within 5% of specification throughout their service life
- B. A significant proportion of in-use curing lights produce irradiance below recommended minimums
- C. LED degradation does not affect irradiance output in devices less than 2 years old
- D. Curing tip contamination has no measurable effect on delivered irradiance

Section 5: Ergonomics and Device Design (Questions 21–24)

21. The Demi Pro curing light weighs approximately:

- A. 75 grams
- B. 115 grams
- C. 185 grams
- D. 240 grams

22. The primary ergonomic benefit of the Demi Pro's 360° rotating curing tip is:

- A. Elimination of the need for light shields during curing
- B. Reduction of heat transfer from the tip to the operator's glove
- C. Allowing optimal light angulation without forced awkward wrist positioning
- D. Increasing irradiance output by optimizing LED alignment

23. Work-related musculoskeletal disorders (MSDs) in dental practitioners most commonly affect which anatomical regions?

- A. Lower back, hip, and knee
- B. Shoulder, neck, and wrist
- C. Elbow, thumb, and ankle
- D. Fingers, forearm, and lower back

24. The Demi Pro offers how many distinct curing modes?

- A. 2 (Standard and High)
- B. 3 (Standard, High, and Ortho)
- C. 4 (Standard, High, Ortho, and Pulse)
- D. 5 (Standard, High, Ortho, Pulse, and Tack)

Section 6: Transillumination (Questions 25–27)

25. Fiber-optic transillumination (FOTI) is classified as which type of diagnostic modality?

- A. Ionizing radiation-based
- B. Thermographic
- C. Noninvasive, radiation-free
- D. Fluorescence-based

26. When using transillumination for caries detection, areas of demineralization appear as:

- A. Brighter zones due to increased light transmission
- B. Unchanged areas with uniform light transmission
- C. Shadows or areas of reduced light transmission
- D. Red-colored zones due to wavelength absorption

27. Transillumination is particularly valuable for identifying which of the following diagnostically challenging conditions?

- A. Class I occlusal caries on primary molars
- B. Vertical root fractures
- C. Open margins on existing PFM crowns
- D. Furcation involvement in multirooted teeth

Section 7: Clinical Protocol Best Practices (Questions 28–30)

28. For conventional composite resins, the maximum increment depth recommended for reliable light-curing is:

- A. 1 mm
- B. 2 mm
- C. 4 mm
- D. 6 mm

29. When should a clinician verify their curing light output using the integrated radiometer?

- A. Only when purchasing a new device
- B. Once per month as part of routine equipment maintenance
- C. At the beginning of each clinical day and after any period of storage, impact, or tip replacement
- D. Exclusively when a restoration fails postoperatively

30. Which of the following BEST summarizes the clinical significance of the Demi Pro's combination of collimated output, PLS technology, and 360° rotating tip?

- A. It reduces the need for incremental composite layering technique
- B. It enables reliable curing performance at clinical distances, thermal safety for the patient, and ergonomic access to all restorative surfaces
- C. It eliminates the requirement for preoperative radiographic assessment
- D. It allows composite restorations to be completed without a rubber dam or matrix band system

CE ANSWER FORM (Email address required for processing)

Name:	Title:	Specialty:
Address:	NPI No.:	
City:	State:	Zip:
AGD Identification No.:		
Email:		
Telephone:	License Renewal Date:	

AGD Code: 250

EDUCATIONAL OBJECTIVES

- Describe the fundamental principles of photopolymerization and the role of light irradiance in achieving adequate composite cure depth and hardness.
- Identify the clinical implications of inadequate light-curing, including marginal failure, secondary caries, and restoration longevity concerns.
- Explain how collimated beam output and curing tip design affect irradiance delivery at clinical distances.
- Evaluate the safety benefits of Periodic Level Shift (PLS) technology in preventing overheating and protecting pulpal and soft tissue integrity.
- Discuss ergonomic design considerations in curing light selection and their relationship to clinician fatigue and procedural efficiency.
- Apply knowledge of multimode curing protocols (Standard, High, Ortho) to appropriate clinical scenarios.
- Recognize the diagnostic applications of transillumination in early-carries detection and fracture assessment.

COURSE EVALUATION

Please evaluate this course using a scale of 5 to 1, where 5 is excellent and 1 is poor.

- | | |
|--|-----------|
| 1. Clarity of objectives | ⑤ ④ ③ ② ① |
| 2. Usefulness of content | ⑤ ④ ③ ② ① |
| 3. Benefit to your clinical practice | ⑤ ④ ③ ② ① |
| 4. Usefulness of the references | ⑤ ④ ③ ② ① |
| 5. Quality of written presentation | ⑤ ④ ③ ② ① |
| 6. Quality of illustrations | ⑤ ④ ③ ② ① |
| 7. Clarity of quiz questions | ⑤ ④ ③ ② ① |
| 8. Relevance of quiz questions | ⑤ ④ ③ ② ① |
| 9. Rate your overall satisfaction with this course | ⑤ ④ ③ ② ① |
| 10. Did this lesson achieve its educational objectives? ☐ Yes ☐ No | |
| 11. Are there any other topics you would like to see presented in the future? _____ | |
| 12. Overall administration of the program | ⑤ ④ ③ ② ① |

COURSE SUBMISSION:

1. Read the entire course.
2. Complete this entire answer sheet in either pen or pencil.
3. Mark only one answer for each question.
4. Fax answer form to 732-303-0555.

For immediate results:

1. Read the entire course.
2. Go to <https://www.dentallearning.net/article/light-curing-outcomes>
3. Log in to your account or register to create an account.
4. Complete course and submit for grading to receive your CE verification certificate.

A minimum score of 70% is required to receive CE credits.

QUIZ ANSWERS

Fill in the circle of the appropriate answer that corresponds to the question on previous page.

- | | |
|-------------------------|-------------------------|
| 1. (A) (B) (C) (D) (E) | 16. (A) (B) (C) (D) (E) |
| 2. (A) (B) (C) (D) (E) | 17. (A) (B) (C) (D) (E) |
| 3. (A) (B) (C) (D) (E) | 18. (A) (B) (C) (D) (E) |
| 4. (A) (B) (C) (D) (E) | 19. (A) (B) (C) (D) (E) |
| 5. (A) (B) (C) (D) (E) | 20. (A) (B) (C) (D) (E) |
| 6. (A) (B) (C) (D) (E) | 21. (A) (B) (C) (D) (E) |
| 7. (A) (B) (C) (D) (E) | 22. (A) (B) (C) (D) (E) |
| 8. (A) (B) (C) (D) (E) | 23. (A) (B) (C) (D) (E) |
| 9. (A) (B) (C) (D) (E) | 24. (A) (B) (C) (D) (E) |
| 10. (A) (B) (C) (D) (E) | 25. (A) (B) (C) (D) (E) |
| 11. (A) (B) (C) (D) (E) | 26. (A) (B) (C) (D) (E) |
| 12. (A) (B) (C) (D) (E) | 27. (A) (B) (C) (D) (E) |
| 13. (A) (B) (C) (D) (E) | 28. (A) (B) (C) (D) (E) |
| 14. (A) (B) (C) (D) (E) | 29. (A) (B) (C) (D) (E) |
| 15. (A) (B) (C) (D) (E) | 30. (A) (B) (C) (D) (E) |

Save time and the environment
by taking this course online.

**If you have any questions,
please email *Dental Learning* at
support@dentallearning.net**

PLEASE PHOTOCOPY ANSWER SHEET FOR ADDITIONAL PARTICIPANTS.

Please direct all questions pertaining to Dental Learning, LLC, or the administration of this course to support@dentallearning.net. **COURSE EVALUATION and PARTICIPANT FEEDBACK:** We encourage participant feedback pertaining to all courses. Please be sure to complete the evaluation included with the course. **INSTRUCTIONS:** All questions have only one answer. Participants will receive confirmation of passing by receipt of a verification certificate. Verification certificates will be processed within two weeks after submitting a completed examination. **EDUCATIONAL DISCLAIMER:** The content in this course is derived from current information and research-based evidence. Any opinions of efficacy or perceived value of any products mentioned in this course and expressed herein are those of the author(s) of the course and do not necessarily reflect those of Dental Learning. Completing a single continuing education course does not provide enough information to make the participant an expert in the field related to the course topic. It is a combination of many educational courses and clinical experience that allows the participant to develop skills and expertise. **COURSE CREDITS:** All participants scoring at least 70% on the examination will receive a CE verification certificate. Dental Learning, LLC, is designated as an Approved PACE Program Provider by the Academy of General Dentistry. The formal continuing education programs of this program provider are accepted by AGD for Fellowship, Mastership, and membership maintenance credit. Please contact Dental Learning, LLC, for current terms of acceptance. Participants are urged to contact their state dental boards for continuing education requirements. Dental Learning, LLC, is a California Provider. The California Provider number is RP5062. **RECORD KEEPING:** Dental Learning, LLC, maintains records of your successful completion of any exam. Please contact our offices for a copy of your continuing education credits report. This report, which will list all credits earned to date, will be generated and mailed to you within five business days of request. Dental Learning, LLC, maintains verification records for a minimum of seven years. Go green, go online to <https://www.dentallearning.net/article/light-curing-outcomes>. © 2026